

Endorsing CONTOURA Vision, Topography-Guided LASIK

Sponsored by Alcon



Endorsing CONTOURA Vision, Topography-Guided LASIK

Four surgeons explain why they use CONTOURA Vision and advise peers on the platform's benefits.



The following article highlights key questions addressed during a roundtable discussion of the role of personalized CONTOURA Vision topography-guided LASIK among today's selection of refractive treatments. Led by Kerry D. Solomon, MD, the panel included Daniel S. Durrie, MD; Ronald R. Krueger, MD; and R. Doyle Stulting, MD, PhD. Video of the roundtable can be viewed at eyetube.net/collections/contoura-pearls.

When surgeons who are performing wavefront-optimized, wavefront-guided, or standard procedures hear about CONTOURA Vision topography-guided LASIK, they want to know how it will help them with routine patients. Why should they consider choosing this approach over what they are already doing?

Ronald R. Krueger, MD: This is the newest form of customized laser vision correction, and it has advantages because the centration is different than what we have done before. Treatment is centered on the apex of the cornea. It is able to regularize the cornea better than any of the procedures available to date, and that gives us the best outcomes for our patients.

I think we have learned a lot from wavefront customization, which has been in many respects the standard way to improve vision for the last 10 years. Now we are taking it an increment further by delivering corneal wavefront. It is the same kind of information we might have about refraction and aberrations, but now we are treating the structure to get the ideal anterior surface of the cornea. If we optimize the corneal surface, we get the kind of results that will improve the optics of the rest of the visual system.

Daniel S. Durrie, MD: If a surgeon who is happy with the results of those procedures asked me, "Why do I need this?" I would ask about the success of their current treatments from a clinical and a business perspective. How have the outcomes of those procedures driven practice growth over the last couple of years? I do not think

we can be satisfied with what we have if our market is not growing. Doctors need next-generation technologies. Yes, they require a financial investment, but our patients deserve a surgeon who keeps driving outcomes up to the next level. CONTOURA Vision topography-guided LASIK allows us to grow this market largely with the great tools we already have.

Clinically—even if the surgeon is living in the high-contrast acuity world, feels satisfied with the 20/20 rates, and does not have a lot of enhancements—I think he or she is missing those people who could go to the next level. We can give patients higher quality night vision, for example, if we use all the tools at our disposal. CONTOURA Vision is another tool that helps us give patients the benefits of increased **quality vision and correction of aberrations** on the cornea that are asymmetry-based.

R. Doyle Stulting, MD, PhD: The scientific or the technical advantage of the CONTOURA Vision topography-guided procedure is that it is more reproducible. It reports more points in the cornea than wavefront does. It allows us to get readings for a larger area of the cornea. It avoids reading from inside the eye. It is immune to abnormalities that are created by lenticular aberrations. Finally, it is not affected by accommodation or affected by various types or intraocular opacities. Those are the reasons that the procedure is theoretically more accurate; results of the clinical trial show that its accuracy produces better visual outcomes.¹

The clinical trial showed that CONTOURA Vision topography-guided LASIK resolves aberrations and visual complaints and has a lower rate of LASIK-induced symptoms compared to other approaches. Is that a key consideration for surgeons considering this approach?

Dr. Stulting: For laser users with a premium practice who want to offer their patients the best outcomes, the results of the CONTOURA Vision trial are illuminating. The trial showed that patients' visual aberrations and complaints were lessened by the treatment.¹ For the first time, we have an opportunity to say that patients who had this treatment have better uncorrected vision

CONTOURA Vision is a trademark of Novartis. ©2017 Novartis. All other brand/product names are trademarks of their respective owners.

©2017 Novartis 08/17 US-WVL-17-E-1879

than their best-corrected vision with glasses or contact lenses before surgery.¹ I want to offer this new option to patients in my practice.

I have been doing LASIK for nearly 20 years. We used to tell patients that the price they would pay for freedom from glasses would be some new visual aberrations. They could expect to have difficulty driving at night and possibly other issues. Now I can legitimately tell patients that based on the results of the clinical trial, it is likely that their visual aberrations will be lessened after treatment.¹ Instead of losing night vision, they may actually gain it. This is something that we have never been able to offer to patients before. I think that is why we need CONTOURA Vision topography-guided LASIK, even for normal eyes.

Kerry D. Solomon, MD: I think surgeons need CONTOURA Vision topography-guided LASIK to give their patients the opportunity for better uncorrected vision. The uncorrected vision results and other visual outcomes data are some of the best prospective study results I have seen to date in the literature. Adding in the reduction of visual symptoms—night vision and others—means we are also making patients' quality of vision better. CONTOURA Vision topography-guided treatment offers better visual outcomes, reduced symptoms, and better quality vision, and therefore it should fit into everyone's armamentarium, whether they choose to use it for all routine patients or select patients who would benefit the most. ■

1. Stulting RD, Fant BS, T-CAT Study Group. Results of topography-guided laser in situ keratomileusis custom ablation treatment with a refractive excimer laser. *J Cataract Refract Surg.* 2016;42(1):11-18.

Kerry D. Solomon, MD

- medical director, Carolina Eyecare Research Institute, Mt. Pleasant, South Carolina, and adjunct clinical professor of ophthalmology at the Medical University of South Carolina in Charleston
- kerry.solomon@carolinaeyecare.com
- financial disclosure: consultant to Alcon

Daniel S. Durrie, MD

- founder of Durrie Vision; clinical professor, director of refractive surgery services and director of the ophthalmic refractive surgery fellowship program at the University of Kansas Medical Center
- ddurrie@DurrieVision.com
- financial disclosure: consultant to Alcon

Ronald R. Krueger, MD

- medical director of refractive surgery at the Cleveland Clinic Cole Eye Institute and professor of ophthalmology at the Cleveland Clinic Lerner College of Medicine of Case Western Reserve University
- krueger@ccf.org
- financial disclosure: consultant to Alcon

R. Doyle Stulting, MD, PhD

- co-founder of the Stulting Research Center, a corneal disease and LASIK specialist, and director of corneal disease and research at Woolfson Eye Institute
- dstulting@woolfsoneye.com
- financial disclosure: consultant to Alcon

WaveLight® Excimer Laser Systems Important Product Information

This information pertains to all WaveLight® Excimer Laser Systems, including the WaveLight® ALLEGRETTO WAVE®, the ALLEGRETTO WAVE® Eye-Q, and the WaveLight® EX500.

Caution: Federal (U.S.) law restricts the WaveLight® Excimer Laser Systems to sale by or on the order of a physician. Only practitioners who are experienced in the medical management and surgical treatment of the cornea, who have been trained in laser refractive surgery (including laser calibration and operation) should use a WaveLight® Excimer Laser System.

Indications: FDA has approved the WaveLight® Excimer Laser systems for use in laser-assisted in situ keratomileusis (LASIK) treatments for:

- the reduction or elimination of myopia of up to - 12.00 D and up to 6.00 D of astigmatism at the spectacle plane;
- the reduction or elimination of hyperopia up to + 6.00 D with and without astigmatic refractive errors up to 5.00 D at the spectacle plane, with a maximum manifest refraction spherical equivalent of + 6.00 D;
- the reduction or elimination of naturally occurring mixed astigmatism of up to 6.00 D at the spectacle plane; and
- the wavefront-guided reduction or elimination of myopia of up to -7.00 D and up to 3.00 D of astigmatism at the spectacle plane.

In addition, FDA has approved the WaveLight® ALLEGRETTO WAVE® Eye-Q Excimer Laser System, when used with the WaveLight® ALLEGRO Topolyzer® and topography-guided treatment planning software for topography-guided LASIK treatments for the reduction or elimination of up to -9.00 D of myopia, or for the reduction or elimination of myopia with astigmatism, with up to -8.00 D of myopia and up to 3.00 D of astigmatism.

The WaveLight® Excimer Laser Systems are only indicated for use in patients who are 18 years of age or older (21 years of age or older for mixed astigmatism) with documentation of a stable manifest refraction defined as ≤ 0.50 D of preoperative spherical equivalent shift over one year prior to surgery, exclusive of changes due to unmasking latent hyperopia.

Contraindications: The WaveLight® Excimer Laser Systems are contraindicated for use with patients who:

- are pregnant or nursing;

- have a diagnosed collagen vascular, autoimmune or immunodeficiency disease;
- have been diagnosed keratoconus or if there are any clinical pictures suggestive of keratoconus;
- are taking isotretinoin (Accutane®) and/or amiodarone hydrochloride (Cordarone®);
- have severe dry eye;
- have corneas too thin for LASIK;
- have recurrent corneal erosion;
- have advanced glaucoma; or
- have uncontrolled diabetes.

Warnings: The WaveLight® Excimer Laser Systems are not recommended for use with patients who have:

- systemic diseases likely to affect wound healing, such as connective tissue disease, insulin dependent diabetes, severe atopic disease or an immunocompromised status;
- a history of Herpes simplex or Herpes zoster keratitis;
- significant dry eye that is unresponsive to treatment;
- severe allergies;
- a history of glaucoma;
- an unreliable preoperative wavefront examination that precludes wavefront-guided treatment; or
- a poor quality preoperative topography map that precludes topography-guided LASIK treatment.

The wavefront-guided LASIK procedure requires accurate and reliable data from the wavefront examination. Every step of every wavefront measurement that may be used as the basis for a wavefront-guided LASIK procedure must be validated by the user. Inaccurate or unreliable data from the wavefront examination will lead to an inaccurate treatment.

Topography-guided LASIK requires preoperative topography maps of sufficient quality to use for planning a topography-guided LASIK treatment. Poor quality topography maps may affect the accuracy of the topography-guided LASIK treatment and may result in poor vision after topography-guided LASIK.

Precautions: The safety and effectiveness of the WaveLight® Excimer Laser Systems have not been established for patients with:

Endorsing CONTOURA Vision, Topography-Guided LASIK

- progressive myopia, hyperopia, astigmatism and/or mixed astigmatism, ocular disease, previous corneal or intraocular surgery, or trauma in the ablation zone;
- corneal abnormalities including, but not limited to, scars, irregular astigmatism and corneal warpage;
- residual corneal thickness after ablation of less than 250 microns due to the increased risk for corneal ectasia;
- pupil size below 7.0 mm after mydriatics where applied for wavefront-guided ablation planning;
- history of glaucoma or ocular hypertension of > 23 mmHg;
- taking the medications sumatriptan succinate (Imitrex®);
- corneal, lens and/or vitreous opacities including, but not limited to cataract;
- iris problems including, but not limited to, coloboma and previous iris surgery compromising proper eye tracking; or
- taking medications likely to affect wound healing including (but not limited to) antimetabolites.

In addition, safety and effectiveness of the WaveLight® Excimer Laser Systems have not been established for:

- treatments with an optical zone < 6.0 mm or > 6.5 mm in diameter, or an ablation zone > 9.0 mm in diameter; or
- wavefront-guided treatment targets different from emmetropia (plano) in which the wavefront calculated defocus (spherical term) has been adjusted;

In the WaveLight® Excimer Laser System clinical studies, there were few subjects with cylinder amounts > 4 D and ≤ 6 D. Not all complications, adverse events, and levels of effectiveness may have been determined for this population.

Pupil sizes should be evaluated under mesopic illumination conditions. Effects of treatment on vision under poor illumination cannot be predicted prior to surgery.

Adverse Events and Complications

Myopia: In the myopia clinical study, 0.2% (2/876) of the eyes had a lost, misplaced, or misaligned flap reported at the 1 month examination.

The following complications were reported 6 months after LASIK: 0.9% (7/818) had ghosting or double images in the operative eye; 0.1% (1/818) of the eyes had a corneal epithelial defect.

Hyperopia: In the hyperopia clinical study, 0.4% (1/276) of the eyes had a retinal detachment or retinal vascular accident reported at the 3 month examination.

The following complications were reported 6 months after LASIK: 0.8% (2/262) of the eyes had a corneal epithelial defect and 0.8% (2/262) had any epithelium in the interface.

Mixed Astigmatism: In the mixed astigmatism clinical study, two adverse events were reported. The first event involved a patient who postoperatively was subject to blunt trauma to the treatment eye 6 days after surgery. The patient was found to have an intact globe with no rupture, inflammation or any dislodgement of the flap. UCVA was decreased due to this event. The second event involved the treatment of an incorrect axis of astigmatism. The axis was treated at 60 degrees instead of 160 degrees.

The following complications were reported 6 months after LASIK: 1.8% (2/111) of the eyes had ghosting or double images in the operative eye.

Wavefront-Guided Myopia: The wavefront-guided myopia clinical study included 374 eyes treated; 188 with wavefront-guided LASIK (Study Cohort) and 186 with Wavefront Optimized® LASIK (Control Cohort). No adverse events occurred during the postoperative period of the wavefront-guided LASIK procedures. In the Control Cohort, one subject undergoing traditional LASIK had the axis of astigmatism programmed as 115 degrees instead of the actual 155 degree axis. This led to cylinder in the left eye.

The following complications were reported 6 months after wavefront-guided LASIK in the Study Cohort: 1.2% (2/166) of the eyes had a corneal epithelial defect; 1.2% (2/166) had foreign body sensation; and 0.6% (1/166) had pain. No complications were reported in the Control Cohort.

Topography-Guided Myopia: There were six adverse events reported in the topography-guided myopia study. Four of the eyes experienced transient or temporary decreases in vision prior to the final 12 month follow-up visit, all of which were resolved by the final follow-up visit. One subject suffered from decreased vision in the treated eye, following blunt force trauma 4 days after surgery. One subject experienced retinal detachment, which was concluded to be unrelated to the surgical procedure.

Clinical Data

Myopia: The myopia clinical study included 901 eyes treated, of which 813 of 866 eligible eyes were followed for 12 months. Accountability at 3 months was 93.8%, at 6 months was 91.9%, and at 12 months was 93.9%. Of the 782 eyes that were eligible for the uncorrected visual acuity (UCVA) analysis of effectiveness at the 6-month stability time point, 98.3% were corrected to 20/40 or better, and 87.7% were corrected to 20/20 or better. Subjects who responded to a patient satisfaction questionnaire before and after LASIK reported the following visual symptoms at a “moderate” or “severe” level at least 1% higher at 3 months post-treatment than at baseline:

visual fluctuations (28.6% vs. 12.8% at baseline).

Long term risks of LASIK for myopia with and without astigmatism have not been studied beyond 12 months.

Hyperopia: The hyperopia clinical study included 290 eyes treated, of which 100 of 290 eligible eyes were followed for 12 months. Accountability at 3 months was 95.2%, at 6 months was 93.9%, and at 12 months was 69.9%. Of the 212 eyes that were eligible for the UCVA analysis of effectiveness at the 6-month stability time point, 95.3% were corrected to 20/40 or better, and 69.4% were corrected to 20/20 or better. Subjects who responded to a patient satisfaction questionnaire before and after LASIK reported the following visual symptoms as “much worse” at 6 months post-treatment: halos (6.4%); visual fluctuations (6.1%); light sensitivity (4.9%); night driving glare (4.2%); and glare from bright lights (3.0%).

Long term risks of LASIK for hyperopia with and without astigmatism have not been studied beyond 12 months.

Mixed Astigmatism: The mixed astigmatism clinical study included 162 eyes treated, of which 111 were eligible to be followed for 6 months. Accountability at 1 month was 99.4%, at 3 months was 96.0%, and at 6 months was 100.0%. Of the 142 eyes that were eligible for the UCVA analysis of effectiveness at the 6-month stability time point, 97.3% achieved acuity of 20/40 or better, and 69.4% achieved acuity of 20/20 or better. Subjects who responded to a patient satisfaction questionnaire before and after LASIK reported the following visual symptoms at a “moderate” or “severe” level at least 1% higher at 3 months post-treatment than at baseline: sensitivity to light (52.9% vs. 43.3% at baseline); visual fluctuations (43.0% vs. 32.1% at baseline); and halos (42.3% vs. 37.0% at baseline).

Long term risks of LASIK for mixed astigmatism have not been studied beyond 6 months.

Wavefront-Guided Myopia: The wavefront-guided myopia clinical study included 374 eyes treated; 188 with wavefront-guided LASIK (Study Cohort) and 186 with Wavefront Optimized® LASIK (Control Cohort). 166 of the Study Cohort and 166 of the Control Cohort were eligible to be followed at 6 months. In the Study Cohort, accountability at 1 month was 96.8%, at 3 months was 96.8%, and at 6 months was 93.3%. In the Control Cohort, accountability at 1 month was 94.6%, at 3 months was 94.6%, and at 6 months was 92.2%.

Of the 166 eyes in the Study Cohort that were eligible for the UCVA analysis of effectiveness at the 6-month stability time point, 99.4% were corrected to 20/40 or better, and 93.4% were corrected to 20/20 or better. Of the 166 eyes in the Control Cohort eligible for the UCVA analysis of effectiveness at the 6-month stability time point, 99.4% were corrected to 20/40 or better, and 92.8% were corrected to 20/20.

In the Study Cohort, subjects who responded to a patient satisfaction questionnaire before and after LASIK reported the following visual symptoms at a “moderate” or “severe” level at least 1% higher at 3 months post-treatment than at baseline: light sensitivity (47.8% vs. 37.2% at baseline) and visual fluctuations (20.0% vs. 13.8% at baseline). In the Control Cohort, the following visual symptoms were reported at a “moderate” or “severe” level at least 1% higher at 3 months post-treatment than at baseline: halos (45.4% vs. 36.6% at baseline) and visual fluctuations (21.9% vs. 18.3% at baseline).

Long term risks of wavefront-guided LASIK for myopia with and without astigmatism have not been studied beyond 6 months.

Topography-Guided Myopia: The topography-guided myopia clinical study included 249 eyes treated, of which 230 eyes were followed for 12 months. Accountability at 3 months was 99.2%, at 6 months was 98.0%, and at 12 months was 92.4%. Of the 247 eyes that were eligible for the UCVA analysis at the 3-month stability time point, 99.2% were corrected to 20/40 or better, and 92.7% were corrected to 20/20 or better. Subjects who responded to a patient satisfaction questionnaire before and after LASIK reported the following visual symptoms as “marked” or “severe” at an incidence greater than 5% at 1 month after surgery: dryness (7% vs. 4% at baseline) and light sensitivity (7% vs. 5% at baseline). Visual symptoms continued to improve with time, and none of the visual symptoms were rated as being “marked” or “severe” with an incidence of at least 5% at 3 months or later after surgery.

Long term risks of topography-guided LASIK for myopia with and without astigmatism have not been studied beyond 12 months.

Information for Patients: Prior to undergoing LASIK surgery with a WaveLight® Excimer Laser System, prospective patients must receive a copy of the relevant Patient Information Booklet, and must be informed of the alternatives for correcting their vision, including (but not limited to) eyeglasses, contact lenses, photorefractive keratectomy, and other refractive surgeries.

Attention: Please refer to a current WaveLight® Excimer Laser System Procedure Manual for a complete listing of the indications, complications, warnings, precautions, and side effects.

* Trademarks are property of their respective owners.