

Hydrodissection Cannula for Laser Cataract Surgery

PRICE	N/A
COMPANY	Sterimedix
PHONE	+44 1527 501480
WEB	www.sterimedix.com
KEY FEATURES	
<ul style="list-style-type: none"> • Enhances hydrodissection after a laser capsulotomy • Angled and flattened on the horizontal plane • Two sideports to improve the flow of fluid 	

Designed in conjunction with a number of leading ophthalmic surgeons in Europe and the United States, this cannula for laser cataract surgery enhances hydrodissection after a laser capsulotomy, according to Sterimedix. The 25-gauge cannula is angled and flattened on the horizontal plane for easier insertion under the capsule. The instrument has two sideports to improve the flow of fluid after dissection to overcome any persistent adherence of the capsule and cortex. Each box is packaged with 10 sterile cannulas.



Richer Glare Disability Test

PRICE	N/A
COMPANY	M&S Technologies
PHONE	(877) 225-6101
WEB	www.mstech-eyes.com
KEY FEATURES	
<ul style="list-style-type: none"> • Available on the Smart System 20/20 • Simulates the effect of oncoming headlights or stepping out into bright sunshine • Validates patients' complaints of poor vision in glare conditions 	

The Richer Glare Disability Test (M&S Technologies) is now available on the Smart System 20/20 via the M&S Glare Testing System. According to the company, the test emits bright glare and then briefly presents a visual acuity letter or symbol, simulating the effect of sudden oncoming headlights or stepping out into bright sunshine. Eye care professionals can use this fast and easy test to validate patients' complaints of poor vision in glare conditions. ■

